CERTEST ASTROVIRUS
One step Astrovirus Card Test

INTENDED USE

The CerTest Astrovirus is a one step coloured chromatographic immunoassay for the qualitative detection of Astrovirus in stool samples.

INTRODUCTION

Astrovirus is a cause of infectious gastroenteritis in infants and young children, also observed in adults. It is transmitted by fecal-oral contact. The main symptoms of viral gastroenteritis are watery diarrhoea and vomiting. The affected person may also have headache, fever, and abdominal cramps ("stomach ache"). In general, the symptoms begin 1 to 2 days following infection with Astrovirus that causes gastroenteritis and may last for 3 days.

PRINCIPLE OF THE TEST

The CerTest Astrovirus is a qualitative immunochromatographic assay for the determination of Astrovirus in feces samples. The membrane is pre-coated with mouse monoclonal antibodies, on the test band region, against viral antigens. During testing, the sample is allowed to react with the coloured conjugate (anti-astrovirus mouse monoclonal antibodies-red microspheres) which was pre-dried on the test. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the coloured particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the coloured conjugate. The mixture continues to move across the membrane to the immobilized antibody placed in the control band region, a GREEN coloured band always appears. The presence of this GREEN band serves as 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) as an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS


- For professional *in vitro* diagnostic use only.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The tests should be discarded in a proper biohazard container after testing.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers and the assay should be done right after collection. The samples can be stored in the refrigerator (2-4°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed, and brought to room temperature before testing.

Specimen preparation (see illustration):

1. Unscrew the cap and use the stick to pick up quite a lot of sample (200-300mg), if the stool sample was liquid take 200-300 µL using a pipette, and add the sample into the stool collection tube.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion.

MATERIALS PROVIDED

- CerTest Astrovirus card-device test
- Instructions for use
- Stool collection tubes with sample diluent


MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer


TEST PROCEDURE

Allow the tests, stool samples and controls to reach to room temperature (15-30°C) prior to testing. Do not open pouches until ready to perform the assay.

1. Proceed to shake the stool collection tube in order to assure good sample dispersion. Cut the end of the cap (3).
2. Remove the CerTest Astrovirus device from its sealed bag just before using.
3. Use a separate stool collection tube and device for each sample or control. Dispense 4 drops or 100 µL into the circular window marked with an arrow (4).
4. Read the result at **10 minutes** (the coloured bands appear).

INTERPRETATION OF RESULTS (please refer to the illustration below)

NEGATIVE: Only one GREEN band appears across the central window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control band, a RED band (test line) also appears in the site marked with the letter T (result line).

INVALID: A total absence of the control coloured band (GREEN) regardless the appearance or not of the result line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the result line region (T) will vary depending on the concentration of antigens in the specimen. However, neither the quantitative value, nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test. A green line appearing in the control region (C) is an internal control. It confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours of opening the sealed bag.
2. An excess of stool sample could cause wrong results (brown bands appear).
3. After one week of infection, the number of viruses in feces is decreasing, making the sample less reactive. Stool samples should be collected within one week of the onset of symptoms.
4. This test provides a presumptive diagnosis for Astrovirus infections. A confirmed infection diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated.

EXPECTED VALUES

Negative results are expected in healthy infants and young children, also in healthy adults.

PERFORMANCE

The evaluation was conducted comparing the results obtained using the CerTest Astrovirus to a commercial available Astrovirus ELISA assay.

Sensitivity

The detection of Astrovirus showed >94% of concordance in sensitivity.

Specificity

The detection of Astrovirus showed a >99% of concordance in specificity.

The use of mouse monoclonal antibodies in the elaboration of CerTest Astrovirus assures high degree of specificity for the detection of this virus.

REFERENCES

1. CUKOR G., and BLACKLOW N. R., "Human Viral Gastroenteritis", *Microbiological Reviews*, Vol. 48 No 2, June 1984, pp. 157-179
2. NEEL K. KRISHNA, B.A., "Identification of Structural Domains Involved in Astrovirus Capsid Biology", *Viral Immunol.* 2005 ; 18(1): 17-26.
3. BON, F. et al. "Prevalence of group A rotavirus, human calicivirus, astrovirus type 40 and 41 infections among children with acute gastroenteritis in Dijon, France." *J. Clin. Microbiol.* 37 No 9 3055-3058 (1999).

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	<i>In vitro</i> diagnostic device		Lot code
	Consult instructions for use		Catalogue number
	Keep dry		Contains sufficient for <n> tests
	Temperature limitation		Manufacturer
	Use by		Sample diluent

CerTest C/ María de Luna 11, nave 16
BIOTEC S.L. E-50018 Zaragoza (SPAIN)
www.cerTEST.es

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