**INTENDED USE**

The CerTest Calprotectin-Lactoferrin is a one step coloured chromatographic immunoassay for the qualitative determination of human calprotectin and lactoferrin in stool samples that reflect intestinal inflammation caused by several pathologies (inflammatory bowel disease, colorectal cancer and some enteropathies). CerTest Calprotectin-Lactoferrin offers a highly sensitive, simple and non-invasive assay to estimate intestinal inflammatory activity, to monitor treatment response and to predict risk of relapse.

**INTRODUCTION**

**Calprotectin** is a neutrophil cytosolic protein with antimicrobial properties, which is present at increased concentration in stool during bowel inflammation. The stability of the protein to degradation keeps it stable in faeces for up to 7 days at room temperature, making it an ideal analyte. Calprotectin is released by activation of leukocytes, giving increased levels in plasma, cerebral spinal fluid, synovial fluid, urine or stools as a consequence of disease in the relevant organ(s). Calprotectin inhibits zinc-dependent enzyme systems, as a result kills microbes and induces apoptosis in normal and cancer cells. In the presence of calcium, calprotectin is a remarkably resistant to proteolytic degradation, and so is stable in stools kept at room temperature for 7 days.

**Lactoferrin** is a glycoprotein component of neutrophil secondary granules, a primary component of the acute inflammatory response, and is released from fecal leukocytes. This protein is resistant to proteolysis in the feces and may serve as a marker of inflammation in the intestine. The major cause of fecal neutrophils in patients with chronic diarrhoea is chronic inflammatory bowel disease of the colon (i.e., Crohn's Disease and ulcerative colitis). Lactoferrin has been also studied as a predictor of infection with invasive enteropathogens in children with diarrhoea. Bacterial inflammatory diarrhea may be caused by *Shigella*, *Salmonella*, *Campylobacter* and *Clostridium difficile*.

These markers can differentiate between functional (for example: irritable bowel syndrome) and organic disease (for example inflammatory bowel disease).

CerTest Calprotectin-Lactoferrin offers a highly sensitive, simple and non-invasive assay to estimate inflammatory disease, to monitor treatment response in these patients and to predict risk of relapse.

**PRINCIPLE OF THE TEST**

The CerTest Calprotectin-Lactoferrin is a qualitative immunochromatographic assay for the determination of human calprotectin and lactoferrin in stool samples. The membrane is pre-coated with antibodies on the test band (result region), against calprotectin and lactoferrin.

During testing, the sample is allowed to react with the coloured conjugates (anti-human calprotectin antibodies-red microspheres and anti-human lactoferrin antibodies) pre-dried on the test. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the coloured particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the coloured conjugate. One or two coloured band will be visible, depend on the calprotectin/lactoferrin content of the sample. These bands are used to interpret the result.

The mixture continues to move across the membrane to the immobilized antibody placed in the control band region; this green coloured band always appears. The presence of this green band serves as 1) verification that sufficient volume is added, 2) proper flow obtained and 3) an internal control for the reagents.

**STORAGE AND STABILITY**

Store as packaged in the sealed pouch at 2-30°C. The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

**PRECAUTIONS**


- For professional *in vitro* diagnostic use only.
- Do not use after expiration date.
- All the specimens should be considered as potentially hazardous and handled as if they were infectious agents.
- The tests should be discarded in a proper biohazard container after testing.

**SPECIMEN COLLECTION AND PREPARATION**

Stool samples should be collected in clean containers and the assay should be done right after collection. The samples can be stored at room temperature or in the refrigerator (2-4°C) for 7 days prior to testing. For longer storage, maximum 6 months, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed, and brought to room temperature before testing.

**Specimen preparation (see illustration):**

- (1) Unscrew the cap and use the stick by introducing four times into the fecal specimen to pick up some sample and adding the sample (10-20mg) into the stool collection tube.
- (2) Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion.

**MATERIALS PROVIDED**

- CerTest Calprotectin-Lactoferrin device tests
- Instructions for use
- Stool collection tubes with sample diluent


**MATERIALS REQUIRED BUT NO PROVIDED**

- Specimen collection container
- Disposable gloves
- Timer


**TEST PROCEDURE**

**Allow the test, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until ready to perform the assay.**

1. Proceed to shake the stool collection tube in order to assure good sample dispersion. Cut the end of the cap (3).
2. Remove the CerTest Calprotectin-Lactoferrin card device from its sealed bag just before using.
3. Use a separate stool collection tube and device for each sample or control. Dispense 4 drops or 100 µL into the circular window marked with an arrow (4).
4. Read the result at **10 minutes** (the coloured bands appear). Do not read the test result after more than 10 minutes.

**INTERPRETATION OF RESULTS (Please refer to the illustration below)**

**NEGATIVE:** Only GREEN band appears across the central window in the site marked with the letter C (control line). That would mean neither active gastrointestinal inflammation is present in the patient. If symptoms continue that would mean irritable bowel syndrome (functional disease).

**Calprotectin (FCP) and Lactoferrin (FLF) POSITIVE:** In addition to the GREEN control band in the control region, a RED band and a BLUE band in the result region) could appear at the same time during the test performance due to a simultaneous presence of elevated fecal calprotectin and lactoferrin in the sample. That would mean intestinal inflammation (organic disease) is present in the patient (caused by ulcerative colitis (UC), Crohn's disease (CD), some type of carcinomas or an enteropathy caused by *Shigella*, *Salmonella*, *Campylobacter* or *Clostridium difficile*).

**Lactoferrin POSITIVE:** In addition to the GREEN control band, a BLUE band (lactoferrin test line) also appears in the site marked with the letter T (result region). That probably would mean gastrointestinal inflammatory pathology cause by CD or UC, and with persistent diarrhoea it probably would be caused by viral or bacterial enteropathies.

**Calprotectin POSITIVE:** In addition to the GREEN control band, a RED band (calprotectin test line) also appears in the site marked with the letter T (result region). That probably would mean gastrointestinal inflammatory pathology cause by CD or UC or some carcinoma.

**INVALID:** A total absence of the control coloured band (GREEN) regardless the appearance or not of the result line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are likely the reasons for control line failure. Review the procedure and repeat the tests using a new test. If the problem persists, discontinue using the test kit and contact you local distributor.

**NOTES ON THE INTERPRETATION OF RESULTS**

The intensity of the red or blue coloured band in the result region (T) will vary depending on the concentration of calprotectin or lactoferrin in the specimen. However, neither the quantitative value, nor the rate of increase in calprotectin or lactoferrin can be determined by this qualitative test.

**QUALITY CONTROL**

Internal procedural controls are included in the test. A green line appearing in the control region (C) is an internal control. It confirms sufficient specimen volume and correct procedural technique.

**LIMITATIONS**

1. The test must be carried out within 2 hours of opening the sealed bag.
2. Stool from patients with non-steroidal anti-inflammatory drug (NSAID) treatment could show positive calprotectin result.
3. Neonatal fecal calprotectin levels have been reported higher than normal children with a median of 167µg/g (range 22-860).
4. An excess of stool sample could result in wrong results (brown bands appear or absence of the control coloured band).
5. Stool from patients with active inflammatory bowel diseases that usually involve significant neutrophilic inflammation of the intestine, such as Crohn's disease and ulcerative colitis, would be positive for fecal calprotectin (the level of calprotectin in Crohn's disease is higher than in ulcerative colitis) and would be positive for fecal lactoferrin. CerTest Calprotectin-Lactoferrin could be sensitive for this diagnosis in patients with chronic diarrhea.
6. Positive results confirm the presence of human calprotectin or/and human lactoferrin in fecal samples; nevertheless, it can be also due to several causes (Inflammatory bowel disease, colorectal cancer and some enteropathies). A positive result should be followed up with additional diagnostic invasive procedures, a colonoscopy and a biopsy in order to confirm the diagnosis and to establish the inflammation extent.
7. Negative results do not exclude inflammation, some diseases such as celiac sprue and microscopic colitis polyps that involve mainly monocuclear inflammation.
8. Lactoferrin is a component of breast milk; the test will be positive in breast fed children and should not be used to evaluate neonates receiving breast milk.

**EXPECTED VALUES**

Some studies established equal or higher 50µg hFCP/g faeces and 10µg hFLF/g faeces as cut-off values to allow detect adult patients with gastrointestinal inflammatory pathology that will require to diagnosis additional diagnostic invasive procedures.

**PERFORMANCE**

**Sensitivity**

A sample containing calprotectin at concentration equal to or higher than 50µg/g faeces and a concentration of lactoferrin equal to or higher than 10µg/g faeces produces positive results when using CerTest Calprotectin-Lactoferrin test.

The detection of human calprotectin with CerTest Calprotectin-Lactoferrin test showed >94% of sensitivity correlation compared to another commercial immunoassay (Calprest® Eurospital). And the detection of human fecal lactoferrin with CerTest Calprotectin-Lactoferrin test showed >99% of sensitivity correlation compared to another commercial immunoassay (IBD EZ VUE®TechLab®).

**Specificity**

The detection of human calprotectin with CerTest Calprotectin-Lactoferrin test showed 93% of specificity correlation compared to Calprest® Eurospital and the detection of human fecal lactoferrin with CerTest Calprotectin-Lactoferrin test showed >99% of specificity correlation compared to IBD EZ VUE®TechLab®.

The CerTest Calprotectin-Lactoferrin test is specific for human fecal calprotectin and human fecal lactoferrin, showing no cross-reaction with bovine lactoferrin.

**REFERENCES**

1. Angriman I. et al. Enzymes in feces: Useful markers of chronic inflammatory bowel disease. *Clinica Chimica Acta* 381 Feb 2007, p. 63-68
2. Quail, M.A. et al. Fecal Calprotectin Complements Routine Laboratory Investigations in Diagnosing Childhood Inflammatory Bowel Disease. *Inflamm Bowel Dis*, Vol 15 No 5; May 2009, p. 756-759
3. Gaya D.R., et al. Faecal calprotectin in the assessment of Crohn's disease activity. *Q J Med* 2005, Vol 98, May 2005, p. 435-441.
4. Langhorst, M.D. et al. Noninvasive Markers in the Assessment of Intestinal Inflammation in Inflammatory Bowel Diseases: Performance of Fecal Lactoferrin, Calprotectin and PMN-Elastase, CRP, and Clinical Indices. *Am. J. Gastroenterol.* 2008; Vol 103, p. 162-169.
5. Gisbert, J.P. et al. Role of Biological Markers in Inflammatory Bowel Disease. *Gastroenterol. Hepatol.* May 2007, 30(3), p. 117-29.

**SYMBOLS FOR IVD COMPONENTS AND REAGENTS**

	<i>In vitro</i> diagnostic device		Lot code
	Consult instructions for use		Catalogue number
	Keep dry		Contains sufficient for <n> tests
	Temperature limitation		Manufacturer
	Use by		Sample diluent

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1. Angriman I. et al. Enzymes in feces: Useful markers of chronic inflammatory bowel disease. *Clinica Chimica Acta* 381 Feb 2007, p. 63-68
2. Quail, M.A. et al. Fecal Calprotectin Complements Routine Laboratory Investigations in Diagnosing Childhood Inflammatory Bowel Disease. *Inflamm Bowel Dis*, Vol 15 No 5; May 2009, p. 756-759
3. Gaya D.R., et al. Faecal calprotectin in the assessment of Crohn's disease activity. *Q J Med* 2005, Vol 98, May 2005, p. 435-441.
4. Langhorst, M.D. et al. Noninvasive Markers in the Assessment of Intestinal Inflammation in Inflammatory Bowel Diseases: Performance of Fecal Lactoferrin, Calprotectin and PMN-Elastase, CRP, and Clinical Indices. *Am. J. Gastroenterol*. 2008; Vol 103, p. 162-169.

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