CERTEST *E.coli* O157:H7
One Step *E.coli* O157:H7 Card test

INTENDED USE

The CerTest *E.coli* O157:H7 is a one step coloured chromatographic immunoassay for the qualitative detection of *Escherichia coli* O157 in stool samples.

INTRODUCTION

Infection with *Escherichia coli* O157:H7 (Enterohemorrhagic *Escherichia coli*) presents with a wide spectrum of clinical manifestations, including asymptomatic carriage, nonbloody diarrhea, hemorrhagic colitis, the hemolytic-uremic syndrome, and thrombotic thrombocytopenic purpura. Not only is *E. coli* O157:H7 an important agent for hemorrhagic colitis, it is also one of the leading causes of bacterial diarrhea. Transmission of *E. coli* O157:H7 is primarily food-borne. Undercooked meat is the most common culprit, dairy products and secondary person-to-person spread are also important. The organism produces at least two Shiga-like toxins. These toxins are thought to have direct pathogenic significance in *E. coli* O157:H7 infection. This infection is usually diagnosed from a positive stool culture, from the presence of Shiga toxins, or both. Timely collection (within 7 days of illness onset) of a stool sample for culture is imperative for a high recovery rate.

CerTest *E. coli* O157:H7 provides a rapid detection of *E. coli* O157:H7 directly from the stool samples.

PRINCIPLE OF THE TEST

The CerTest *E. coli* O157:H7 is a qualitative immunochromatographic assay for the determination of *E. coli* O157:H7 in stool samples. The membrane is pre-coated with mouse polyclonal antibodies, on the test band region, to recognize this serotype (polyclonal antibodies against O157:H7 antigens).

During testing, the sample is allowed to react with the coloured conjugate (anti-*E. coli* O157:H7 antibodies-red microspheres) which was pre-dried on the test. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the coloured particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the coloured conjugate (result region). The mixture continues to move across the membrane to the immobilized antibody placed in the control band region, a GREEN coloured band always appears. The presence of this GREEN band serves as 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) as an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS


- For professional *in vitro* diagnostic use only.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The tests should be discarded in a proper biohazard container after testing.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers and the assay should be done right after collection. The samples can be stored in the refrigerator (2-4 °C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed, and brought to room temperature before testing.

Specimen preparation (see illustration):

- (1) Unscrew the cap and use the stick to pick up a little sample (approx. 100mg), if the stool sample was liquid take 100 µL using a pipette, and add the sample into the stool collection tube.
- (2) Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion.

Note: Because it is possible for the *E. coli* content in the stool to be too small, a negative result cannot be considered to be the same as a negative EHEC finding. It is recommended to carry out enrichment in TSB medium with Mitomycin C additive overnight on a shaker at 37°C.

MATERIALS PROVIDED

- Device tests (CerTest *E. coli* O157:H7)
- Instructions for use
- Stool collection tubes with sample diluent

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer


TEST PROCEDURE

Allow the tests, stool or enrichment samples and controls to reach to room temperature (15-30°C) prior to testing. Do not open pouches until ready to perform the assay.


1. Proceed to shake the stool collection tube in order to assure good sample dispersion. Cut the end of the cap (3).
2. Remove the CerTest *E. coli* O157:H7 device from its sealed bag just before using.
3. Use a separate stool collection tube and device for each sample or control. Dispense exactly 4 drops or 100 µL into the circular window marked with an arrow (4).

In case the test did not run due to solid particles fallen into the round window, stir the sample added or dispense a drop of extraction buffer until seeing the liquid running through the reaction zone.

- 4.- Read the result at **10 minutes** (the coloured bands appear).

INTERPRETATION OF RESULTS (please refer to the illustration below)

NEGATIVE: Only one GREEN band appears across the central window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control band, a distinguishable RED band also appears in the site marked with the letter T (result line).

INVALID: A total absence of the control coloured band (GREEN) regardless of the appearance or not of the result line (RED). Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the tests with a new test. If the problem persists, discontinue using the test kit and contact you local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the result line region (T) will vary depending on the concentration of antigen present in the specimen. However, neither the quantitative value, nor the rate of increase in antigen can be determined by this qualitative test.

QUALITY CONTROL

A green line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours of opening the sealed bag.
2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
3. Some stool samples can decrease the intensity of the control green line.
4. A negative result is not meaningful because it is possible the *E.coli* content in the stool sample to be too small. A *E. coli* determination should be carried out on a sample from a enrichment culture.
5. This test provides a presumptive diagnosis of *E. coli* infections. A confirmed infection diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated must be based in the correlation of the results with further clinical observations.

PERFORMANCE

CerTest *E. coli* O157:H7 was evaluated to determine: sensibility in stool samples directly and in selective enrichment cultures, specificity with producers organisms of Shiga toxins, non-Shiga toxins producers and other *Enterobacteriaceae* species.

Sensitivity

The detection of *E. coli* O157:H7 showed a 100% of concordance in sensitivity.

Specificity

The detection of *E. coli* O157:H7 showed a 85% of concordance in specificity.

PPV showed a 70% and NPV showed a 100%.

Exactitude

The detection of *E. coli* O157:H7 showed a 94% of concordance in exactitude.

REFERENCES

1. THOMPSON, J., HODGE, D. and BORCZYK, A.; "Rapid Biochemical Test to Identify Verocytotoxin-Positive Strains of *Escherichia coli* Serotype O157"; *Journal of Clinical Microbiology*, Oct. 1990, Vol. 28, No. 10, pp 2165-2168
2. VALLANCE B.A. and FINLAY B.B., "Exploitation of host cells by enteropathogenic *Escherichia coli*", *PNAS*, August 2000, Vol. 97, No. 16, pp. 8799-8806
3. BLANCO, M. et al. " *Escherichia coli* Verotoxigenicos (ECVT) en España: ECVT O157:H7 y NO-O157 en humanos y alimentos. El ganado bovino y ovino como reservorio. Técnicas para detección de ECVT" Laboratorio de Referencia de *Escherichia coli* (LREC), Departamento de Microbiología y Parasitología, Facultad de Veterinaria, Universidad de Santiago de Compostela, 27002 Lugo

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	In vitro diagnostic device		Batch code
	Consult instructions for use		Catalogue number
	Keep dry		Contains sufficient for <n> tests
	Temperature limitation		Manufacturer
	Use by		Sample diluent

CerTest BIOTEC S.L. C/María de Luna 11, nave 16 E-50018 Zaragoza (SPAIN) www.cer-test.es

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