CERTEST ENTEROVIRUS

One step Enteroviruses Card Test in stool samples

INTENDED USE

The CerTest Enterovirus is a one step coloured chromatographic immunoassay for the qualitative detection of Enteroviruses (VP1 peptide) in stool samples. Enteroviruses may cause a wide variety of pathological symptoms and enteroviral infections that affect especially young children. CerTest Enterovirus allows for early detection of enterovirus, that is very interesting in epidemiology.

INTRODUCTION

Enteroviruses belong to the *Picornaviridae* family and are the smallest, non-enveloped viruses known to infect both humans and animals. The spread of enteroviral infections is mainly by the faecal-oral and oral-oral route, but also through direct contact with secretions from ophthalmic and dermal lesions. Water, food and soil contaminated by infected faeces are an exogenous infection source which creates many opportunities for the transfer of the infection, and cause an epidemic outbreak in a short period of time. The average incubation period for enteroviral contagious is from 3–10 to 30 days.

The virus, after replicating and breaking the gastrointestinal tract barrier, is transmitted via the blood stream to every organ of the body. Enteroviruses reveal tropism towards organs like the heart, skin, and in particular the central nervous system. It has been proved that infected people excrete large quantities of the virus in faeces for a period of even 16 weeks.

Traditional division organizes this taxonomic group into the subgenera polioviruses, coxsackieviruses (group A, B), echoviruses and a group of enteroviruses marked according to their serotype number (66–71 and newly identified 73–75, 77, 78). In 2003, the International Committee on Taxonomy of Viruses created a new taxonomy classification. Enteroviruses henceforth were divided into 5 groups of species based on their molecular properties (HEV-A, HEV-B, HEV-C, HEV-D).

PRINCIPLE OF THE TEST

The CerTest Enterovirus is a qualitative immunochromatographic assay for the determination of Enteroviruses in fecal samples. The membrane is pre-coated with mouse monoclonal antibodies, on the test band region, against viral antigens (VP1 peptide). This antibody reacts with an epitope on the VP1 peptide which is highly conserved within the enterovirus group. During testing, the sample is allowed to react with the coloured conjugate (anti-VP1 peptide mouse monoclonal antibodies-red microspheres) which was pre-dried on the test. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the coloured particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the coloured conjugate. The mixture continues to move across the membrane to the immobilized antibody placed in the control band region, a GREEN coloured band always appears. The presence of this GREEN band serves as 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) as an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS


- For professional *in vitro* diagnostic use only.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The tests should be discarded in a proper biohazard container after testing.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers and the assay should be done right after collection. The samples can be stored in the refrigerator (2-4°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed, and brought to room temperature before testing.

Specimen preparation (see illustration):

- (1) Unscrew the cap and use the stick to pick up a little sample, if the stool sample was liquid take 100 µL using a pipette, and add the sample into the stool collection tube.
- (2) Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion.

MATERIALS PROVIDED

- Device tests (CerTest Enterovirus)
- Instructions for use
- Stool collection tubes with sample diluent


MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer


TEST PROCEDURE

Allow the tests, stool samples and controls to reach to room temperature (15-30°C) prior to testing. Do not open pouches until ready to perform the assay.

1. Proceed to shake the stool collection tube in order to assure good sample dispersion. Cut the end of the cap (3).
2. Remove the CerTest Enterovirus device from its sealed bag just before using.
3. Use a separate stool collection tube and device for each sample or control. Dispense 4 drops or 100 µL into the circular window marked with an arrow (4).
- 4.- Read the result at **10 minutes** (the coloured bands appear).

INTERPRETATION OF RESULTS (please refer to the illustration below)

NEGATIVE: Only one GREEN band appears across the central window in the site marked with the letter C (control line). A negative result indicates no Enterovirus infection. A confirmed diagnosis should be made after all clinical and laboratory findings.

POSITIVE: In addition to the GREEN control band, a RED band (test line) also appears in the site marked with the letter T (result line). This result indicates a recent and contagious infection with Enterovirus.

INVALID: A total absence of the control coloured band (GREEN) regardless the appearance or not of the result line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the result line region (T) will vary depending on the concentration of antigens in the specimen. However, neither the quantitative value, nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test. A green line appearing in the control region (C) is an internal control. It confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours of opening the sealed bag.
2. An excess of stool sample could cause wrong results (brown bands appear).
3. After one month of infection, the number of viruses in feces is decreasing, making the sample less reactive. Stool samples could be collected previously to the onset of symptoms or also at 24-48 hours.
4. This test provides a presumptive diagnosis for *Enterovirus* group infections. A confirmed infection diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated and should be determine the serotype of enterovirus to established the type of disease.
5. If the patient has been recently vaccinated (for example against poliovirus), it would be appeared a positive result.

EXPECTED VALUES

Negative results are expected in healthy infants and young children, also in healthy adults. Infections caused by the enteroviruses are most likely to occur during the summer and fall.

PERFORMANCE

Sensitivity and specificity

It was studied some stool samples from patients with suspicious of enterovirus infection. All samples were evaluated by culture. It was performed an evaluation using CerTest Enterovirus and the results showed >99% of sensitivity and >99% of specificity. The antibodies used to elaborate this test recognise epitopes present in antigens found in stool patients.

Cross-reactivity

It was performed an evaluation to determine the cross reactivity of CerTest *Enterovirus*. There is not cross reactivity with other pathogen intestinal virus occasionally present in faeces:

- Rotavirus	- Astrovirus
- Adenovirus	- Hepatitis A

REFERENCES

1. RAJTAR, B. et al. "Enteroviruses in water environment-a potential threat to public health". *Ann. Agric. Environ. Med.*, 2008, No 15, p. 199-203.
2. BAEK, K. et al. "Molecular and Epidemiological Characterization of Enteroviruses Isolated in Chungnam, Korea from 2005-2006." *J. Microbiol. Biotechnol.*, 2009, 19(9), p. 1055-1064.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	<i>In vitro</i> diagnostic device		Lot code
	Consult instructions for use		Catalogue number
	Keep dry		Contains sufficient for <n> tests
	Temperature limitation		Manufacturer
	Use by		Sample diluent

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