CERTEST *H. pylori*
Rapid One step *Helicobacter pylori* Card test

INTENDED USE

The CerTest *H. pylori* is a one step coloured chromatographic immunoassay for the qualitative detection of *Helicobacter pylori* in faeces.

INTRODUCTION

Helicobacter pylori (*H. pylori*) is a spiral-shaped bacterium that is found in the gastric mucous layer or adherent to the epithelial lining of the stomach. *H. pylori* causes more than 90% of duodenal ulcers and up to 80% of gastric ulcers.

The importance of *Helicobacter pylori* testing has increased greatly since the strong correlation between the presence of bacteria and confirmed gastrointestinal diseases (stomach and duodenum) like gastritis, peptic ulcer disease and gastric carcinoma.

PRINCIPLE OF THE TEST

The CerTest *H. pylori* is a qualitative immunochromatographic assay for the determination of *Helicobacter pylori* in faeces samples. The membrane is pre-coated with monoclonal antibodies, on the test band region, against *H. pylori* antigens.

During testing, the sample is allowed to react with the coloured conjugate (anti-*H. pylori* monoclonal antibodies-red polystyrene micro spheres) which was pre-dried on the test strip. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the coloured particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the coloured conjugate. The mixture continues to move across the membrane to the immobilized antibody placed in the control band region, a green coloured band always appears. The presence of this green band serves as 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) as an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS


- For professional *in vitro* diagnostic use only.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The tests should be discarded in a proper biohazard container after testing.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers and the assay should be done right after collection. The samples can be stored in the refrigerator (2-4 °C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed, and brought to room temperature before testing.

Specimen preparation (see illustration):

- (1) Take out the top of the stool collection tube and use the stick to pick up enough quantity of sample, introduce the stick 4 times and collect approx. 250 mg. Add the sample into the stool collection tube.
- (2) Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion.

Note: if the stool sample was liquid take 250 µL using a pipette.

MATERIALS PROVIDED

- Device tests (CerTest *H. pylori* card test)
- Instructions for use
- Stool collection tubes with sample diluent

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer


TEST PROCEDURE

Allow the tests, stool samples and controls to reach to room temperature (15-30°C) prior to testing. Do not open pouches until ready to perform the assay.


1. Proceed to shake the stool collection tube in order to assure good sample dispersion. Cut the end of the top (3).
2. Remove the CerTest *H. pylori* device from its sealed bag just before using.
3. Use a separate stool collection tube and device for each sample or control. Dispense exactly 4 drops or 100 µL into the circular window marked with an arrow, avoiding to add solid particles with the liquid (4).

In case the tests did not run due to solid particles fallen into the round window, stir the sample added or dispense a drop of extraction buffer until seeing the liquid running through the reaction zone.

- 4.- Read the result at **10 minutes** (the coloured bands appear).

INTERPRETATION OF RESULTS (please refer to the illustration below)

NEGATIVE: Only one GREEN band appears across the central window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control band, a RED band (test line) also appears in the site marked with the letter T (result line).

INVALID: A total absence of the control coloured band regardless the appearance or not of the result line. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the result line region (T) will vary depending on the concentration of antigens in the specimen. However, neither the quantitative value, nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test. A green line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique. A clear background is an internal negative background control. If the test is working properly, the background in the result area should be clear and not interfere with the ability to read the result.

LIMITATIONS

1. The test must be carried out within 2 hours of opening the sealed bag.
2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
3. Some stool samples can decrease the intensity of the control line.
4. This test provides a presumptive diagnosis of *Helicobacter pylori* infections. A confirmed infection diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated must be based in the correlation of the results with further clinical observations.

PERFORMANCE

Sensitivity and specificity

It was studied some patients with the same as *H. pylori* infection symptoms. For all patient, it was performed an evaluation using CerTest *H. pylori* and a commercial available ELISA *H. pylori* assay to detect *H. pylori* infection. The results are showed:


	EIA Test	
	Sensitivity	Specificity
CerTest <i>H. pylori</i>	> 94%	> 99%

The use of a mouse monoclonal antibody in CerTest *H. pylori* assures high degree of specificity for the detection of this bacteria. The antibodies used to elaborate the *H. pylori* Test recognise epitopes present in the antigen found in stool of patients, as well as in preparations from the bacteria cultures *in vitro*. Sonicated *Helicobacter pylori* extract from different commercial samples reacts with CerTest *H. pylori*.

REFERENCES

- 1- Bruce E. Dunn, Hartley Cohen & Martin J. Blaser. *Helicobacter pylori*. Clin. Microbiol. Rev. **10** (4), 720-741, Oct. (1997)
- 2- Martin J. Blaser. *Helicobacter pylori and gastric diseases*. BMJ; **316**: 1507-1510 (1998).
- 3- John L. Telford, Antonello Covacci, Rino Rappuoli & Paolo Ghiara. *Immunology of Helicobacter pylori infections*. Current Opinion in Immunology, **9**: 498-503 (1997).

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

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