## BioTracer™ HBsAg Rapid Card

※ Product code : 10112(30T/Kit)/10135(100T/Kit)

### INTENDED USE

**BioTracer™ HBsAg Rapid Card** is intended for the qualitative determination of HBsAg in human serum or plasma.

### EXPLANATION OF THE TEST

**Hepatitis B** is a systemic disease primarily involving the liver. The complex antigen found on the surface of HBV is called HBsAg. The presence of HBsAg in serum or plasma is an indication of an active Hepatitis B infection, either acute or chronic. Hundreds of millions of people, most of them in regions with poor medical care, are chronically infected with the virus and face an elevated risk of acquiring liver cancer. **BioTracer™ HBsAg Rapid Card** use monoclonal anti-HBsAg (from mouse) to capture and monoclonal anti-HBsAg(from mouse)-colloidal gold conjugate to detect the presence of HBsAg in human serum or plasma specimen.

### PRINCIPLE OF THE METHOD

**BioTracer™ HBsAg Rapid Card** is based on the principle of an immunochromatography *in vitro* test for the qualitative determination of HBsAg. When the sample is added to sample pad, it moves to the conjugate pad and resuspends anti-HBsAg antibody-conjugated gold complex that is dried on the conjugate pad. The mixture moves along the membrane by capillary action and reacts with anti-HBsAg antibody that is immobilized on the test reaction zone. If HBsAg is present enough in the sample, a colored band in the test reaction zone is appeared. If there is no HBsAg or not sufficient in the sample, the area will remain colorless. The sample continues to move to the control reaction zone and forms a red or purple color, indicating the test is working properly and the result is valid.

### CONTENTS

1. Test device
2. Package insert

### SPECIMEN COLLECTION AND PREPARATION

1. Centrifuge whole blood to get plasma/serum specimens.
2. Testing should be performed immediately after the specimens have been collected. Do not leave the specimens at room temperature for prolonged period. Specimens may be stored and tested at 2-8 °C for up to three days. For long term storage, specimens should be kept below -20 °C.
3. Specimens containing precipitate may yield inconsistent test results. Such specimens must be clarified prior to assaying.

### TEST PROCEDURE

1. All specimens and test devices should be prepared with warm condition, that is, for 15~30 min at room temperature before testing.
2. Remove the test device from the sealed pouch, and place it on a clean and flat surface.
3. Transfer 100µl of serum or plasma into the sample well(S).
4. Wait for the red or purple band in the result window to appear after the test begins to work.
5. **Interpret test results at 15-20 minutes.**


※ Caution : Perform the test immediately after removing the test device from the foil pouch.

### READING AND INTERPRETATION OF RESULT

1. A color band will appear at the control reaction zone(C) of the result window to show that the test is working properly.
2. A color band of the test reaction zone (T) indicates the test result.

#### NEGATIVE:

The presence of only one band of the control reaction zone(C) indicates a negative result.


#### POSITIVE :

The presence of two color bands (T) band and (C) band within the result window regardless of which band appears first indicates a positive result.

#### INVALID :

If control(C) band is not appeared in the result window after performing test, the result is considered invalid. The directions may not have been followed correctly or the test device may have been deteriorated. It is recommended that the specimen be re-tested with the new test device.

### FOLLOW-UP ACTION

As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

### LIMITATIONS OF THE METHOD

A negative result does not preclude the possibility of infection with HBV. Other clinically available tests are required if questionable results are obtained.

### DETECTION LIMIT

The detection limit of the **BioTracer™ HBsAg Rapid Card** is 1ng/ml.

### PERFORMANCE CHARACTERISTICS

A study was performed using 256 positive and negative serum specimens. Each specimen was assayed with the **BioTracer™ HBsAg Rapid Card** and a commercially available HBsAg test(ELISA) according to the respective package insert.

	ELISA	BioTracer™ HBsAg Rapid Card		Total
		Positive	Negative	
	Positive	163	1	164
	Negative	0	92	92
	Total	163	93	256

The sensitivity is 99.4%(163/164) and the specificity is more than 99.9%(92/92).

### WARNING AND PRECAUTIONS

1. For *in vitro* diagnostic use only.
2. Do not eat or smoke while handling specimens.
3. Wear protective gloves while handling specimens. Wash hands thoroughly afterwards.
4. Avoid splashing or aerosol formation.
5. Clean up spills thoroughly using an appropriate disinfectant.
6. Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials, as if they were infectious waste, in a biohazard container.
7. Do not use the test kit if the packing is damaged or the bottle is broken.

### STORAGE AND SHELF-LIFE

**BioTracer™ HBsAg Rapid Card** should be stored at 2-30 °C (36-86 °F). The test device is sensitive to humidity as well as to heat. Do not use it beyond the expiry date, 18 months from manufacturing.

### LITERATURE REFERENCES

1. Aubuchon JP, Sandler SG, Fang CT, et al, "American Red Cross Experience With Routine Testing for Hepatitis B Core Antibody," *Transfusion*, 1989, 29(3):230-2.
2. Chambers LA and Popovsky MA, "Decrease in Reported Posttransfusion Hepatitis. Contributions of Donor Screening for Alanine Aminotransferase and Antibodies to Hepatitis B Core Antigen and Changes in the General Population," *Arch Intern Med*, 1991, 151(12):2445-8.
3. Bortolotti F, Calzia R, Cadrobbi P, et al, "Long-Term Evolution of Chronic Hepatitis B in Children With Antibody to Hepatitis B Antigen," *J Pediatr*, 1990, 116(4):552-5.
4. Edwards MS, "Hepatitis B Serology – Help in interpretation" *Pediatr Clin North Am*, 1988,35:503-15

Manufactured by  
**NanoEnTek Inc.**  
 851-14, Seohaeh-ro, Paltan-myeon, Hwaseong-si,  
 Gyeonggi-do, 18531, Korea  
 Tel : +82-2-6220-7942 / Fax : +82-2-6220-7999  
 ivdst@nanoentek.com  
 www.nanoentek.com

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