NanoEntek

HCV Rapid Card

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※ Product code : 11012(30T/Kit) / 11035(100T/Kit)

INTENDED PURPOSE

BioTracer™ HCV Rapid Card is intended for the qualitative determination of antibody against HCV in human whole blood, serum or plasma.

EXPLANATION OF THE TEST

Hepatitis C Virus (HCV) is a major cause of acute hepatitis and is widespread among injecting drug users, with up to 90% being positive for antibodies to HCV. HCV infection progresses with an indolent course and remains benign in the first few years, but the majority (>50%) of HCV-infected individuals eventually develop chronic HCV infection. HCV was first reported in 1989. It is now known that HCV is an enveloped virus that contains a positive strand RNA genome of about 9400 nucleotides. Antibody against HCV is found in most patients who are infected with HCV. **BioTracer™ HCV Rapid Card** use recombinant HCV structural protein(core) and non-structural protein(NS3, NS4, NS5) to capture and protein A-colloidal gold conjugate to detect the antibodies against HCV in serum, plasma or whole blood. This test is intended for professional use as an aid in the detection of antibodies against HCV.

PRINCIPLE OF THE METHOD

BioTracer™ HCV Rapid Card is based on the principle of an immunochromatography *in vitro* test for the qualitative determination of antibody against HCV. When the sample is added to sample pad, it moves to the conjugate pad and resuspends protein A conjugated gold complex that is dried on the conjugate pad. The mixture moves along the membrane by capillary action and reacts with HCV antigen that is immobilized the test reaction zone. If antibody against HCV is present enough in the sample, a colored band in the test reaction zone is appeared. If there is no antibody against HCV or not sufficient in the sample, the area will remain colorless. The sample continues to move to the control reaction zone and forms a red or purple color, indicating the test is working properly and the result is valid.

CONTENTS

1. Test device
2. Sample diluent
3. Package insert

SPECIMEN COLLECTION AND PREPARATION

1. Whole blood specimen collection
 - 1) Whole blood is collected in syringe or evacuated tube containing the anticoagulant.
 - 2) Whole blood specimens should be tested immediately after collection. In the case of storing at 2-8°C, it should be tested within 24 hours.
2. Plasma / Serum specimen collection
 - 1) Plasma or serum specimens should be tested immediately after collection.
 - 2) Do not leave the specimens at room temperature for prolonged period. Specimens may be stored at 2-8°C for up to 3 days. For long term storage, specimens should be kept below -20°C.

※ Specimens containing precipitate may yield inconsistent test results. Such specimens must be clarified prior to assaying.

TEST PROCEDURE

Allow the test devices and specimens to room temperature before removal from the pouch and prior to use.

1. Remove the test device from the sealed pouch, and place it on a clean and flat surface.
2. Transfer 10 μ l of serum(or plasma) or whole blood into the sample well(S). And then add 80-100 μ l (4-5 drops) of sample diluent into the sample well(S). On treating the whole blood, one should transfer that sample and sample diluent into the sample well almost simultaneously, otherwise, it may cause clogging of the sample well(S).
3. As the test begin to work, you will see purple (or red) color to move across the result window in the center of the test device.
4. **Interpret test result at 10 minutes.** Do not interpret test result after 15minutes.

※ Caution: Perform the test immediately after removing the test device from the foil pouch.

READING AND INTERPRETATION OF RESULT

1. A color band will appear at the control reaction zone(C) of the result window to show that the test is working properly.
2. A color band of the test reaction zone (T) indicates the test result.

NEGATIVE:

The presence of only one band of the control reaction zone(C) indicates a negative result.


POSITIVE :

The presence of two color bands ("T" band and "C" band) within the result window regardless of which band appears first indicates a positive result.

INVALID :

If control(C) band is not appeared in the result window after performing test, the result is considered invalid. The directions may not have been followed correctly or the test device may have been deteriorated. It is recommended that the specimen be re-tested with the new device.

FOLLOW-UP ACTION

As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

LIMITATIONS OF THE METHOD

A negative result does not preclude the possibility of infection with HCV. Other clinically available tests are required if questionable results are obtained.

PERFORMANCE CHARACTERISTICS

1. Sensitivity and Specificity

A study was performed using 354 positive and negative serum specimens. Each specimen was assayed with the **BioTracer™ HCV Rapid Card** and a commercially available HCV test (ELISA) according to the respective package insert.

The results are summarized in the following tables.

		BioTracer™ HCV Rapid Card		Total
		Positive	Negative	
ELISA	Positive	152	2	154
	Negative	2	198	200
Total		154	200	354

The sensitivity is 98.7%(152/154) and the specificity is 99%(198/200).

2. Reproducibility

- 1) Within run performance test was determined by one analyst with ten devices each of 3 different lots for 4 different positive specimens. There was no variation within the test devices of each lot and between three different lots.
- 2) Between run performance test was determined by three analyst with 3 different lots for 4 different positive specimens. There was no variation between different analyst.

WARNING AND PRECAUTIONS

1. For professional *in vitro* diagnostic use only.
2. Do not eat or smoke while handling specimens.
3. Wear protective gloves while handling specimens. Wash hands thoroughly afterwards.
4. Do not use test kit if the packing is damaged or the seal is broken.
5. Avoid splashing or aerosol formation.
6. Clean up spills thoroughly using an appropriate disinfectant.
7. Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials, as if they were infectious waste, in a biohazard container.

STORAGE AND SHELF-LIFE

BioTracer™ HCV Rapid Card should be stored at 1~30°C (34~86°F). The test device is sensitive to humidity as well as to heat. Do not use it beyond the expiry date, 18 months from manufacturing.

LITERATURE REFERENCES

1. Choo, Q.L., G. Kuo, A.J. Weiner, L.R. Overby, D.W. Bradley, and M. Houghton. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989; 244:359
2. Kuo, G., Q.L. Choo, H.J. Alter, and M. Houghton. An assay for circulating antibodies to a major etiologic Virus of human non-A, non-B hepatitis. *Science* 1989; 244:362
3. Van der Poel, C. L., H.T.M. Cuyper, H.W. Reesink, and P.N.Lelie. Confirmation of hepatitis C Virus infection by new four-antigen recombinant immunoblot assay. *Lancet* 1991; 337:317

Manufactured by
NanoEntek Inc.
 851-14, Seohaeng-ro, Paltan-myeon, Hwaseong-si,
 Gyeonggi-do, 18531, Korea
 Tel : +82-2-6220-7942 / Fax : +82-2-6220-7999
 ivdst@nanoentek.com
 www.nanoentek.com