CERTEST Salmonella typhi+paratyphi

ONE STEP Salmonella
typhi+paratyphi COMBO CARD
TEST

CERTEST BIOTEC S.L.

CERTEST Salmonella typhi+paratyphi

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ENGLISH

INTENDED USE

The CerTest *Salmonella typhi+paratyphi* one step combo card test is a coloured chromatographic immunoassay for the qualitative detection of *Salmonella typhi* and *Salmonella paratyphi* in stool samples to detect a possible infection caused by *Salmonella typhi* and/or *Salmonella paratyphi* in humans.

CerTest *Salmonella typhi+paratyphi* combo card test offers a highly sensitive, simple and non-invasive screening assay to make a presumptive diagnosis of *Salmonella typhi* and/or *Salmonella paratyphi* infection.

INTRODUCTION

Enteric fever (typhoid and paratyphoid fever) is a major human bacterial infection. Although the disease is not common in industrialised countries, it remains an important and persistent health problem in developing nations. That enteric fever is a major public health problem in those countries, with *Salmonella enterica serovar typhi* (*Salmonella typhi*) the most common aetiological agent but with an apparently increasing number of cases due to *Salmonella paratyphi*. Because risk factors such as poor sanitation, lack of a safe drinking water supply and low socio economic conditions in resource-poor countries are amplified by the evolution of multidrug resistant *salmonellae* with reduced susceptibility to fluoroquinolone, which is associated with increased mortality and morbidity.

In Europe, *Salmonella typhi* and *Salmonella paratyphi* infections occurs among travelers returning from disease endemic areas.

Enteric fever caused by *Salmonella paratyphi* is indistinguishable from that caused by *Salmonella typhi*. This fever usually develops one to three weeks after exposure and varies in severity. Symptoms include high fever, weakness, lethargy, muscle pain, headache, loss of appetite and diarrhea or constipation. Pink spots appear on the chest; examinations will usually reveal enlargement of the liver and spleen. In severe cases, symptoms of altered mental status and meningitis (fever, stiff neck, seizures) have been reported.

TEST PRINCIPLE

The CerTest *Salmonella typhi+paratyphi* combo card test is based on the principle of a qualitative immunochromatographic assay for the determination of *Salmonella typhi* and *Salmonella paratyphi* in stool samples.

Strip A consists of: a) membrane pre-coated with antibodies on the test line (T) in the result window, against *Salmonella typhi* antigen and immobilised antibodies on the control line (C) against a specific protein, b) label/sample absorbent pad sprayed with test label solution of mouse monoclonal antibodies anti-*Salmonella typhi* conjugated to red treated polystyrene latex and control label solution of specific binding protein conjugated to green treated polystyrene latex, forming coloured conjugate complexes.

Strip B consists of: a) membrane pre-coated with antibodies on the test line (T) in the result window, against *Salmonella paratyphi* antigen and immobilised antibodies on the control line (C) against a specific protein, b) label/sample absorbent pad sprayed with test label solution of mouse monoclonal antibodies anti-*Salmonella paratyphi* conjugated to red treated polystyrene latex and control label solution of specific binding protein conjugated to green treated polystyrene latex, forming coloured conjugate complexes.

During testing, in each strip, the sample is allowed to react with the coloured conjugates and then the mixture moves upward on the membrane by capillary action. As the sample flows through the test membrane, the coloured particles migrate. In the case of a positive result, the specific antibodies present on the membrane will capture the coloured conjugate. A red coloured band will be visible, in **A** and/or **B** depends on the *Salmonella typhi* or *Salmonella paratyphi* content of the sample. These bands are used to interpret the result.

The mixture continues to move across the membrane to the immobilized antibody placed in the control band region, a green coloured band always appears. The presence of this green band serves as 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) as an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new combo card test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.


SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The assay should be performed immediately after collection. The samples can be stored in the refrigerator (2-8 °C) for 1-2 days prior to testing. For longer storage, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Ensure only the amount needed is thawed because of freezing and defrosting cycles are not recommended. Homogenise stool sample as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

1. Take out the cap of the stool collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once in 4 different parts of the stool sample (2), to collect large-quantity faecal sample (approx. 250mg) and add to the stool collection tube.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

MATERIALS

MATERIALS PROVIDED

- CerTest *Salmonella typhi+paratyphi* combo card test
- Instructions for use
- Stool collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED


- Specimen collection container
- Disposable gloves
- Timer

TEST PROCEDURE

Allow tests, stool samples and controls to reach to room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the stool collection tube in order to assure good sample dispersion. Cut the end of the cap (4)
2. Remove the CerTest *Salmonella typhi+paratyphi* combo card test from its sealed bag just before using.
3. Use a combo card test for each sample or control. Dispense 4 drops or 100µL in the circular window marked with the letter **A** (5) and 4 drops or 100µL, using the same tube, in the circular window marked with the letter **B** (6). Avoid adding solid particles with the liquid.
4. Read the results at **10 minutes** (the coloured bands appear). Do not read the test result after more than 10 minutes.

In case the tests did not run due to solid particles fallen into the round window, stir the sample added or dispense a drop of extraction buffer until seeing the liquid running through the reaction zone.

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


Salmonella paratyphi Strip B-procedure

Add 4 drops in B-circular window

INTERPRETATION OF THE RESULTS (please refer to the illustration below)

- A: Green → Negative *Salmonella typhi*
B: Green → Negative *Salmonella paratyphi*
- A: Green/Red → Positive *Salmonella typhi*
B: Green/Red → Positive *Salmonella paratyphi*
- A: Green/Red → Positive *Salmonella typhi*
B: Green → Negative *Salmonella paratyphi*
- A: Green → Negative *Salmonella typhi*
B: Green/Red → Positive *Salmonella paratyphi*

INVALID RESULTS

-

	A (<i>Salmonella typhi</i>)	B (<i>Salmonella paratyphi</i>)	Interpretation of the results
1.	- GREEN	- GREEN	No infection caused by <i>Salmonella typhi</i> and <i>Salmonella paratyphi</i> .
2.	+ GREEN-RED	+ GREEN-RED	Infection caused by <i>Salmonella typhi</i> and <i>Salmonella paratyphi</i> .
3.	+ GREEN-RED	- GREEN	Infection caused by <i>Salmonella typhi</i> .
4.	- GREEN	+ GREEN-RED	Infection caused by <i>Salmonella paratyphi</i> .
5.	Any other result	Any other result	Invalid result. Invalid result either A or B, we recommend repeating the assay using the same sample with another combo card test.

INVALID: total absence of the control coloured band (GREEN) regardless the appearance or not of the results lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the test using a new combo test. If the problem still persists, discontinue using the test kit and contact your local distributor.

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NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results windows will vary depending on the concentration of antigen of *Salmonella typhi* and *Salmonella paratyphi* present in the specimen. However, neither the quantitative value nor the rate of increase in antigen can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the tests. The GREEN line appearing in the control line (C) in the results window is an internal control. It confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours of opening the sealed bag.
2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
3. Some stool samples can decrease the intensity of the control green line.
4. Positive frozen and thawed samples could cause wrong results. Use fresh stool samples if possible.
5. A negative result is not meaningful because it is possible the *Salmonella typhi* and/or the *Salmonella paratyphi* content in the stool sample to be too small. A *Salmonella typhi/Salmonella paratyphi* determination should be carried out on a sample from an enrichment culture.
6. This test provides a presumptive diagnosis of *Salmonella typhi* and *Salmonella paratyphi* infections. A confirmed infection diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated must be based in the correlation of the results with further clinical observations.

EXPECTED VALUES

An estimated 22 million cases of typhoid fever and 200.000 related deaths occur worldwide each year; and additional 6 million cases of paratyphoid fever are estimated to occur annually.

Approximately 400 cases of typhoid fever and 150 cases of paratyphoid fever are reported to CDC each year among persons with onset of illness in the United States, most of whom are recent travellers.

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity.

It was performed an evaluation using a *Salmonella typhi* and *Salmonella paratyphi* cultures. The correlation between assays, Singlepath®Salmonella (Merck) and CerTest *Salmonella typhi+paratyphi* was >99%.

Different bacteria culture extract dilutions were tested directly in the sample diluent or spiked in a negative stool samples. These preparations were evaluated (positive and negative samples) using both assays (Singlepath®Salmonella (Merck) and CerTest *Salmonella typhi+paratyphi*) in accordance with the kit instructions.

The results showed a high sensitivity and specificity using CerTest *Salmonella typhi+paratyphi*.

The antibodies used to elaborate this test recognise epitopes found in stool of patients, as well as in preparations from the bacteria cultures *in vitro*.

<i>Salmonella typhi</i>	IC test: CerTest <i>Salmonella typhi+paratyphi</i>	IC test: Singlepath®Salmonella (Merck)		
		+	-	
	+	10	0	10
	-	0	25	25
		10	25	35

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Cross reactivity

It was performed an evaluation to determine the cross reactivity of CerTest *Salmonella typhi+paratyphi* combo card test. There is not cross reactivity with common intestinal pathogens, other organisms and substances occasionally present in faeces:

<i>Campylobacter jejuni</i>	<i>Helicobacter pylori</i>	<i>Shigella dysenteriae</i>
<i>Campylobacter coli</i>	<i>Listeria monocytogenes</i>	<i>Shigella flexneri</i>
<i>Clostridium difficile</i>	<i>Staphylococcus aureus</i>	<i>Shigella sonnei</i>
<i>E. coli</i> O157:H7	<i>Shigella boydii</i>	<i>Yersinia enterocolitica</i>

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