CERTEST SWINE FLU

One step Influenza A card for the detection of Influenza type A from nasal human specimens and from nasal and tissue or stool pig specimens

INTENDED USE

The CerTest Swine Flu card is a one step coloured chromatographic immunoassay for the qualitative detection of Influenza type A antigens. It can be used directly with nasal swabs or nasal wash or nasal aspirate specimens from human and nasal swabs, lung lavage fluid, lung tissue or stool samples from pigs.

INTRODUCTION

Since swine influenza viruses were first isolated, they have become an economically important cause of respiratory disease in pigs throughout the world, and a human public health risk.

Swine influenza viruses are most commonly of the H1N1 subtype, but other subtypes are also circulating in pigs (e.g., H1N2, H3N1, H3N2). Pigs can also be infected with avian influenza viruses and human seasonal influenza viruses as well as swine influenza viruses. The H3N2 swine virus was thought to have been originally introduced into pigs by humans. Sometimes pigs can be infected with more than one virus type at a time, which can allow the genes from these viruses to mix. This can result in an influenza virus containing genes from a number of sources, called a "reassortant" virus.

Although swine influenza viruses are normally species specific and only infect pigs, they do sometimes cross the species barrier to cause disease in humans. These interspecies infections are most likely to occur when people are in close proximity to pigs, such as in swine production barns, livestock exhibits at fairs, and slaughterhouses.

The clinical signs/symptoms of influenza in pigs and people are remarkably similar, with fever, lethargy, lack of appetite and coughing prominent in both species.

PRINCIPLE OF THE TEST

The CerTest Swine Flu card is a qualitative immunochromatographic assay for the determination of Influenza type A antigens. The membrane is pre-coated with mouse monoclonal antibodies, on the test band region, to recognize these antigens.

During testing, the sample is allowed to react with the coloured conjugate (anti-type A mouse monoclonal antibodies-red microspheres) which was pre-dried on the test. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the coloured particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the coloured particles. Different coloured lines will be visible, depending upon the type of antigens contents of the sample. These lines are used to interpret the result.

The mixture continues to move across the membrane to the immobilized antibody placed in the control band region, a GREEN coloured band always appears. The presence of this GREEN band serves as 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) as an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- CerTest Swine Flu card is intended for professional *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not use test devices if the aluminium pouch has been damaged during storage.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new swab must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.

SPECIMEN COLLECTION AND PREPARATION

Human samples

- Nasal swab specimen: Collect specimen with a sterile swab. Process the swab as soon as possible after collecting the specimen.
- Nasal Wash or Aspirate specimen. Samples should be processed as soon as possible after collection. The samples can be stored in the refrigerator (2-4 °C) for 8 hours prior to testing.

Pig samples

- Nasal swabs: Collect specimen with a sterile swab. Process the swab as soon as possible after collecting the specimen.
- Lung lavage fluid: Saline in leak-proof sterile containers. The fluid can be stored in the refrigerator (2-4°C) for 8 hours prior to testing.
- Lung tissue macerated: Take some affected portions of lung tissue and triturate into a zip-lock plastic bag.
- Stool samples: Stool samples could be collected in clean containers. The assay should be done right after collection. The samples can be stored in the refrigerator (2-4 °C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed, and brought to room temperature before testing.

MATERIALS PROVIDED

- CerTest Swine Flu card tests
- Reagent B (Swine Flu)
- Testing tubes
- Sterile swabs
- Disposable pipettes
- Tissue/Stool collection tube with buffer
- Instructions for use

MATERIALS REQUIRED BUT NO PROVIDED


- Disposable gloves
- Timer

TEST PROCEDURE

Allow the tests, swabs and controls to reach to room temperature (15-30°C) prior to testing. Do not open the package until ready to perform the assay. Only bring to room temperature the number of tests required to assay before opening it.


- Procedure A using nasal swab samples:

1. Add 15 drops or 450µl(1) Reagent B and immediately put the swab into the tube.
2. Mix the solution by rotating the swab forcefully against the side of the tube at least 10 times. Best results are obtained when the specimen is vigorously extracted in the solution (2). Extract as much liquid as possible from the swab, squeezing the sides of the tube as the swab is withdrawn. Discard the swab.
3. Remove the CerTest Swine Flu card device from its sealed bag just before using.
4. Place the test on a flat surface. Use a separate pipette and device for each sample or control.
5. Read the result at 10 minutes.

- Procedure B using nasal wash, aspirate samples and lung lavage pig sample:

1. Add 6 drops (1) of the nasal wash, aspirate samples or lung lavage sample with a pipette and 3 drops (2) of Reagent B in a testing tube.
2. Remove the CerTest Swine Flu card-device from its sealed bag just before using. Place the test on a flat surface.
3. Mix the solution with the pipette at least 10 times. Best results are obtained when the specimen is vigorously extracted in the solution (3).
4. Dispense 4 drops or 100 µL from the testing tube, into the circular window marked with an arrow (4).
5. Read the result at 10 minutes.

- (1) Add 6 drops of the collected sample
- (2) Add 3 drops of Reagent B
- (3) Mix the solution with the pipette at least 10 times

Wait 10 minutes

- Procedure C: using tissues and stool samples from pigs.


1. Take out the top of the tissue/stool collection tube (1).
2. **Tissue sample:** Take approximately 100mg triturate lung tissue and add it into the tissue/stool collection tube (2).

Stool sample: Use the stick to pick up a little sample (approx. 100mg), if the stool sample was liquid take 100 µL using a pipette, and add the sample into the stool collection tube (2).


3. Close the tube with the buffer and sample. Shake the tube in order to assure good sample dispersion (3).
4. Remove the CerTest Swine Flu card-device from its sealed bag just before using.
5. Cut the end of the top (4). Use a separate tissue/stool collection tube and device for each sample or control. Dispense exactly 4 drops or 100 µL into the circular window marked with an arrow, avoiding to add solid particles with the liquid (5).

In case the tests did not run due to solid particles fallen into the round window, stir the sample added or dispense a drop of extraction buffer until seeing the liquid running through the reaction zone.

7. Read the result at 10 minutes (the coloured bands appear)

INTERPRETATION OF RESULTS (please refer to the illustration below)

NEGATIVE: Only one GREEN band appears across the central window in the site marked with the letter C (control line).

Influenza A POSITIVE: In addition to the GREEN control band, a RED band (type A test line) also appears in the site marked with the letter T (results lines).

INVALID: A total absence of the control coloured band (GREEN) regardless of the appearance or not of the result line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the result line region (T) will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value, nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test. A green line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours of opening the sealed pack.
2. This test provides a presumptive diagnosis for Influenza A infections. A confirmed infection diagnosis should only be made by a physician/veterinarian after all clinical and laboratory findings have been evaluated.
3. A negative result may be obtained if the specimen is inadequate or antigen concentration is below the sensitivity of the test. Therefore, it is recommended that all negative CerTest Swine Flu test results undergo confirmatory testing using other method and/or by cell culture.

PERFORMANCE

Sensitivity

Different virus extract dilutions were tested directly in the sample diluent or spiked in a negative specimen in accordance with the kit instructions.

The detection of Influenza type A with CerTest Swine Flu test showed >99% of sensitivity compared with another commercial rapid test.

Specificity

The use of mouse monoclonal antibodies in the elaboration of CerTest Swine Flu test assures high degree of specificity for the detection of Influenza type A antigens.

The detection of Influenza type A with CerTest Swine Flu test showed >99% of specificity compared with another commercial rapid test.

REFERENCES

1. Daniel E. Noyola, et al. Comparison of a New Neuraminidase Detection Assay with an Enzyme Immunoassay, Immunofluorescence, and Culture for Rapid Detection of Influenza A and B Viruses in Nasal Wash Specimens. *Journal of Clinical Microbiology*, Mar. 2000, p. 1161-1165 Vol. 38, No. 3
2. Newman A.P. et al. Human Case of Swine Influenza A (H1N1) Triple Reassortant Virus Infection, Wisconsin. *Emerging Infectious Diseases*, Sept 2008, Vol 14, No. 9, pp1470-1472
3. Donruehthai Sreta et al. Pathogenesis of swine influenza virus (Thai isolates) in weaning pigs: an experimental trial *Virology Journal*, 2009, Vol 6, No 34.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	In vitro diagnostic device		Batch code
	Consult instructions for use		Catalogue number
	Keep dry		Contains sufficient for <n> tests
	Temperature limitation		Manufacturer
	Use by		For IVD performance evaluation only

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1. Daniel E. Noyola, et al. Comparison of a New Neuraminidase Detection Assay with an Enzyme Immunoassay, Immunofluorescence, and Culture for Rapid Detection of Influenza A and B Viruses in Nasal Wash Specimens. Journal of Clinical Microbiology, Mar. 2000, p. 1161-1165 Vol. 38, No. 3.
2. Newman A.P. et al. Human Case of Swine Influenza A (H1N1) Triple Reassortant Virus Infection, Wisconsin. Emerging Infectious Diseases, Sept 2008, Vol 14, No. 9, pp1470-1472
3. Donruehthai Sreta et al. Pathogenesis of swine influenza virus (Thai isolates) in weaning pigs: an experimental trial Virology Journal, 2009, Vol 6, No 34.

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