CERTEST
Shigella
dysenteriae
+ STX

ONE STEP
Shigella dysenteriae +
Shiga toxin (STX) COMBO
CARD TEST

CERTEST BIOTEC S.L.

CERTEST *Shigella dysenteriae*+STX

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ENGLISH

INTENDED USE

CerTest *Shigella dysenteriae*+STX one step combo card test is a coloured chromatographic immunoassay for the simultaneous qualitative detection of *Shigella dysenteriae* type 1 and Shiga toxin (STX) produced by *Shigella dysenteriae* in stool samples and *Shigella dysenteriae* suspected colonies in stool culture.

CerTest *Shigella dysenteriae*+STX combo card test offers a simple and highly sensitive screening assay to make a presumptive diagnosis of *Shigella dysenteriae* infection (shigellosis) and it could be used to identify of suspected isolates of *Shigella dysenteriae* from selective media (stool culture).

INTRODUCTION

The four species of the genus *Shigella*; *S. dysenteriae*, *S. flexneri*, *S. boydii* and *S. sonnei* cause a wide spectrum of illness from watery diarrhoea to fulminant dysentery. The low infectious inoculum, as few as 10 organisms, render *Shigella* highly contagious. Shigellosis therefore occurs as an endemic disease in populations characterized by over-crowding, poor housing, poor sanitation and inadequate water supply. The predominant serogroups of *Shigella* occurring in a region also appears to be related to the socioeconomic development; and evidence also indicates that the severity of shigellosis is related to the infecting serogroup.

S. dysenteriae serotype 1 (also known as *S. dysenteriae* type 1 or *S. dysenteriae* 1), has been recognized as the major cause of epidemic dysentery for nearly 100 years. Pandemics of Shiga dysentery have spread across Central America, Bangladesh, South Asia and Central and East Africa.

All four species of *Shigella*, *Shigella dysenteriae*, *Shigella flexneri*, *Shigella boydii* and *Shigella sonnei*, contain large plasmids, which are required for the invasion of bacteria into epithelial cells, as well as virulence loci encoded on the chromosome. However, only one species of *Shigella*, *S. dysenteriae* serotype 1, has been demonstrated to produce a toxin, called Shiga toxin (STX).

Infections with Shiga toxin-producing bacteria are known to cause bloody diarrhea (hemorrhagic colitis). In some patients, most often children and the elderly, infections are accompanied by a more serious disease, the hemolytic uremic syndrome (HUS), which is characterized by acute renal failure, thrombocytopenia and microangiopathic hemolytic anemia.

TEST PRINCIPLE

CerTest *Shigella dysenteriae*+STX is based on the principle of a qualitative immunochromatographic assay for the determination of *Shigella dysenteriae* and Shiga toxin (STX) produced by *Shigella dysenteriae* in stool samples and *Shigella dysenteriae* suspected colonies in stool culture.

Strip A consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Shigella dysenteriae* and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*Shigella dysenteriae*) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

Strip B consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against STX and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-STX) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is *Shigella dysenteriae* positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-*Shigella dysenteriae* monoclonal antibodies-red polystyrene microspheres) in the strip A, and if the sample is STX positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-STX monoclonal antibodies-red polystyrene microspheres) in strip B, which were previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane, by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-*Shigella dysenteriae* antibodies present on the membrane of strip A (test line) and the anti-STX antibodies present on the membrane of strip B (test line) capture the coloured conjugate and the red line will be visible in both strips. These bands are used to interpret the result.

If the sample is negative, there is no *Shigella dysenteriae* and STX antigens presence and yet, the antigens may be present in a concentration lower than the detection limit values, for which the reaction will not take place with any red-coloured conjugate complex. The anti-*Shigella dysenteriae* and the anti-STX antibodies present on the membranes (test lines) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red lines will not appear.

Whether the sample is positive or not, in both strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on both membranes will capture control green-conjugate

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complex and both control lines will always appear. The presence of these green lines serve as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new combo card test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the CerTest *Shigella dysenteriae*+STX test. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.


SPECIMEN COLLECTION AND PREPARATION

Stool samples

Samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Ensure only the amount needed is thawed because of freezing and defrosting cycles are not recommended. Homogenise stool samples as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):


1. Take out the cap of the stool collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2), to collect faecal sample (approx. 125mg) and add it to the stool collection tube. For liquid samples, add approx. 125µL in the collection tube using a micropipette.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

Suspected *Shigella dysenteriae* colonies in stool culture

The Hektoen Enteric Agar (HEA) is the method of choice for the isolation of *Shigella dysenteriae* (aerobic atmosphere, 24 hours/37°C). After 24 hours of incubation in Hektoen the typical *Shigella dysenteriae* colonies will be green or bluish without any black centre.

1. Examine Hektoen plates after 1 day incubation. Select *Shigella dysenteriae* typical colonies. Take out the cap of the collection tube (1). Use the stick or an inoculating needle to pick up 3 or 4 suspected *Shigella dysenteriae* colonies (2) and add them to the collection tube.
2. Close the tube with the diluent and suspected colonies. Shake the tube in order to assure good sample dispersion (3).

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

Salmonella Shigella Agar (SS) is another method for the isolation of *Shigella dysenteriae* (aerobic atmosphere, 24-48 hours/35-37°C). After the incubation period in SS the typical *Shigella dysenteriae* colonies will be transparent without any black center. The SS agar is a very inhibitor that sometimes prevents the growth of certain species of *Shigella*. We recommend using this medium unless accompanied by another inhibitor (HEA).

1. Examine SS plates after incubation period. Select *Shigella* typical colonies. Take out the cap of the collection tube (1). Use the stick or an inoculating needle to pick up 3 or 4 suspected *Shigella dysenteriae* colonies (2) and add them to the collection tube.
2. Close the tube with the diluent and suspected colonies. Shake the tube in order to assure good sample dispersion (3).

Shigella dysenteriae colonies in Salmonella Shigella Agar

MATERIALS

MATERIALS PROVIDED

- CerTest *Shigella dysenteriae*+STX combo card tests
- Instructions for use
- Stool collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer
- Hektoen Enteric Agar (HEA)/Salmonella Shigella Agar (SS)

TEST PROCEDURE


Allow tests, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the stool collection tube in order to assure good sample dispersion.
2. Remove CerTest *Shigella dysenteriae*+STX combo card test from its sealed bag just before using it.
3. Take the stool collection tube, cut the end of the cap (4) and dispense 4 drops in the circular window marked with the letter A (5) and 4 drops, using the same tube, in the circular window marked with the letter B (6). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.


If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of extraction diluent until seeing the liquid running through the reaction zone.

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

Shigella dysenteriae Strip A-procedure


Shiga toxin Strip B-procedure

INTERPRETATION OF THE RESULTS (please refer to the illustration below)

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	A (<i>Shigella dysenteriae</i>)	B (Shiga toxin (STX))	Interpretation of the results
1.	- GREEN	- GREEN	There is no <i>Shigella dysenteriae</i> presence. No infection caused by <i>Shigella dysenteriae</i> .
2.	+ GREEN-RED	+ GREEN-RED	There is <i>Shigella dysenteriae</i> and Shiga toxin (STX) presence. Infection caused by <i>Shigella dysenteriae</i> .
3.	+ GREEN-RED	- GREEN	There is <i>Shigella dysenteriae</i> presence. Infection caused by <i>Shigella dysenteriae</i> .
4.	- GREEN	+ GREEN-RED	There is Shiga toxin (STX) presence. Infection caused by <i>Shigella dysenteriae</i> or by enterohaemorrhagic <i>E. coli</i> (EHEC).
5.	Any other result	Any other result	Invalid result, either A or B, we recommend repeating the assay using the same sample with another test.

INVALID: Total absence of any control coloured line (GREEN) regardless the appearance or not of the result lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test. The green lines appearing in the control lines (C) in the results windows are internal controls, which confirm sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
3. The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
4. CerTest *Shigella dysenteriae*+STX should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
5. Positive results determine the presence of *Shigella dysenteriae* and/or STX of *Shigella dysenteriae* in faecal samples; nevertheless, it can be due to several causes and/or species. A positive result should be followed up with additional laboratory techniques (biochemical and serological methods) to confirm the results. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
6. A negative result is not meaningful because of it is possible the antigens concentration in the stool samples is lower than the detection limit value. If the symptoms or situation persist, a *Shigella dysenteriae* determination should be carried out, on a sample from an enrichment culture.
7. Shiga toxin (STX) produced by *Shigella dysenteriae* is very similar to one produced by EHEC (verotoxin 1). Verotoxin 1 differs from Shiga toxin by one to seven aminoacids differences. CerTest *Shigella dysenteriae*+STX may give a positive result with Verotoxin1 produced by EHEC. *Shigella dysenteriae* type 1 and EHEC determination should be carried out on selective culture and with biochemical analysis.

EXPECTED VALUES

Bacillary dysentery caused by *Shigella dysenteriae* type 1 is a major public health problem in many developing countries. There are four known species of *Shigella* which are pathogenic and infection with *S. dysenteriae* type 1 usually progresses to the most severe form of dysentery with life-threatening complications. The most common underlying cause of death in fatal shigellosis is severe colitis and the immediate associated causes are septicemia and pneumonia. Infection with *S. dysenteriae* type 1 can occur in an epidemic form, and *Shigella*-contaminated food and drink are often the source of epidemic spread.

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PERFORMANCE CHARACTERISTICS

Analytical sensitivity (detection limit)

Detection limit value of CerTest *Shigella dysenteriae*+STX is 3.125×10^4 CFU/mL for *Shigella dysenteriae* and 20.66ng/mL for STX.

Clinical sensitivity and specificity

An evaluation, with prepared faecal samples, was performed using bacteria cultures *in vitro* to determine the detection of *Shigella dysenteriae*. The results were as follows:

Clinical sensitivity and specificity (Strip A)

IC test: CerTest <i>Shigella dysenteriae</i> +STX (<i>Shigella dysenteriae</i>)	Culture confirmation			Total
		+	-	
+		5	0	5
-		0	10	10
Total		5	10	15

	Sensitivity	Specificity	PPV	NPV
CerTest <i>Shigella dysenteriae</i> +STX vs Culture (<i>Shigella dysenteriae</i>)	>99%	>99%	>99%	>99%

Clinical sensitivity and specificity (Strip B)

IC test: CerTest <i>Shigella dysenteriae</i> +STX (STX)	Culture confirmation			Total
		+	-	
+		3	0	3
-		0	40	40
Total		3	40	43

	Sensitivity	Specificity	PPV	NPV
CerTest <i>Shigella dysenteriae</i> +STX vs Culture (STX)	>99%	>99%	>99%	>99%

The results showed a high sensitivity and specificity to detect *Shigella dysenteriae* and Shiga toxin (STX) using CerTest *Shigella dysenteriae*+STX.

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *Shigella dysenteriae*+STX; no cross reactivity against gastrointestinal pathogens occasionally present in faeces:

<i>Escherichia coli</i> O157:H7	<i>Salmonella typhi</i>	<i>Salmonella enteritidis</i>
<i>Helicobacter pylori</i>	<i>Salmonella typhimurium</i>	<i>Staphylococcus aureus</i>
<i>Listeria monocytogenes</i>	<i>Salmonella paratyphi</i>	<i>Yersinia enterocolitica</i>

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- OPINTAN J.A., et al. "Distribution of Serogroups and Serotypes of Multiple Drug Resistant *Shigella* Isolates" *Ghana Medical Association* 2007.
- Welton I., et al. "Isolation of *Shigellae*. VI. Performance of Media with Stool Specimens. *Applied Microbiology*, Sept. 1968, Vol. 16, No. 9, p. 1387-1393.
- Marie A. McDonough and Joan R. Buterton "Spontaneous tandem amplification and deletion of the Shiga toxin operon in *Shigella dysenteriae* 1. *Molecular Microbiology* (1999) 34 (5), 1058-1069.
- Jonas Bergan, Anne Berit Dyve Lingelem, Roger Simm, Tore Skotland, Kirsten Sandvig "Shiga toxins" *Toxicon* 60 (2012) 1085-1107, Elsevier Ltd.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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