

Biosealer CR 6-PS with Ergonomic Sealing Handle.



Biosealer CR 6-PS with Ergonomic Sealing Handle.

Biosealer CR 6-PS is a battery operated sealing unit for PVC tubes, blood bag tubes and sets for plasmaferes. It uses cadmium free batteries with a capacity for more than 500 seals. Biosealer CR 6-PS is made for mobile use and is easy to carry hanging on your shoulder or in your pocket.

Biosealer CR 6-PS is comprised of four units:

- Coaxial cable.
- Power unit.
- Charger.
- Carrying strap.

Note:

The Sealing Handle fits even to Biosealer CR 4 and Biosealer CR 6.

Technical data:

Batteries	12V/2Ah.
Sealing time	Adjustable 0,5-5 seconds (normal 1,5 seconds).
Frequency	40,68 Mhz crystal controlled.
Output	100W/50 Ohm maximum effect.
Dimensions	167 x 141 x 35 mm incl. RF-outlet.
Weight	1,3 kg incl batteries.

This device complies with requirements in MDD 93/42 EEC (CE marked).

LJUNGBERG & KÖGEL AB
Sweden

www.ljungberg-kogel.com