



# Intact PTH ELISA

## Intact Parathyroid Hormone Assay

Intact PTH ELISA is a direct label 2/site intended for the quantitative determination of Intact PTH in human serum.

The measurement of intact parathyroid hormone is the most specific for the diagnosis of primary hyperparathyroidism

**Diametra is part of IDS Group**

[www.diametra.com](http://www.diametra.com)



# Intact PTH ELISA

## Intact Parathyroid Hormone Assay

### Features and Benefits

- *The test offers a high sensitivity (0,49 pg/mL) and simple 90 + 15 min. streptavidin coated microplate assay procedure.*
- *2 internal controls included*
- *Diametra Intact PTH ELISA was compared to a commercially available kit.*  
 *$Y = 1,0124X + 2,7005$   $r^2 = 0,996$  nr. samples = 110*
- *Sensitivity: 0,49 pg/mL*
- *Normal range: 9 - 94 pg/mL*
- *Analytical range: 0,49 -1000 pg/mL*
- *Precision*
  - *Inter Assay : < 8,2 %*
  - *Intra Assay < 6,3%*

### Ordering information

Cat. Number

DKO 157

Intact PTH ELISA

Diametra Srl - Head Office  
20090 Segrate - Milano - Italy  
Via Calabria 15  
Ph +39.02.21.39.184 - Fax +39.02.21.33.354  
V.A.T. - IT 12585340156  
info\_diametra.@idsplc.com

R&D and Manufacturing  
06034 Spello - Perugia - Italy  
Via Pozzuolo 14  
Ph +39 0742 24 851 - Fax +39 0742 31 61 97

www.diametra.com