


ORGENTEC'S **iVISION** Automated analysis of your immunoblots



YOU WISH TO CARRY OUT AUTOMATED AND STANDARDIZED ANALYSES OF YOUR ORGENTEC IMMUNOBLOTS.

- **iVISION Scanware** ensures objective reading of your immunoblots
- reproducible evaluation of the results
- visual comparison with a control strip is not needed

YOU NEED A RELIABLE RESULT.

- **iVISION Scanware** juxtaposes an image of the test strip with the test result for each antigen: graphically, as an index, and as a semiquantitative gradation
- displays a report summary for up to 16 strips ...
- ... as well as individual, patient-specific reports
- strip-specific calibration of the patient bands with the cut-off band
- use of three control bands ensures the accuracy of the test procedure and result

YOUR LABORATORY BUDGET DEMANDS AN ECONOMICAL, EFFICIENT SOLUTION.

- **iVISION Scanware** may be installed at up to five work stations with a single license
- PC currently in use is easily enabled with a USB dongle
- no specific scanner required (minimum resolution 300dpi); you choose the solution most economical for your laboratory
- from test strip to result within a few minutes: scan the template with the immunoblot strips – read the results
- strip archive: search by date, patient, or strip



ORGENTEC Diagnostika GmbH
Carl-Zeiss-Straße 49-51
55129 Mainz
Germany

Phone: +49 6131/9258-0
Fax: +49 6131/9258-58
orgentec@orgentec.com
www.orgentec.com